# PolyMem 敷料 made easy

Wounds International 2015 年 5 月 www.woundsinternational.com

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1. Department of Health. *Our health, our care, our say: a new direction for community services*, Norwich: Stationery Office (Cm 6737), 2006. Available at: [www.official-documents.gov.uk/document/cm67/6737/6737.pdf](http://www.official-documents.gov.uk/document/cm67/6737/6737.pdf)
2. Genet N, Boerma WG, Kringos DS et al. Home care in Europe: a systematic literature review. *BMC Health Serv Res* 2011; 11: 207
3. Vowden P. Hard to heal wounds Made Easy. *Wounds International*, 2011. <http://bit.ly/1EYUGOE>
4. White R. Hard to heal wounds: results of an international survey. *Wounds UK* 2011; 7(4): 22-8
5. Cutting KC, Weigand C, Vowden P. Wound inflammation and the role of dressings. *Wounds International* 2015; 6(2): 44-46
6. European Wound Management Association. Position document: *Pain and trauma at wound dressing changes*. MEP Ltd: London, 2002
7. Beitz AJ, Newman A, Kahn AR et al. A polymeric membrane dressing with antinociceptive properties: analysis with a rodent model of stab wound secondary hyperalgesia. *J Pain* 2004; 5(1): 38-47
8. Kahn AR, Sessions RW, Apasova EV. A superficial cutaneous dressing inhibits pain, inflammation and swelling in deep tissues. Poster. World Pain Conference, July 2000.
9. Davies SL, White RJ. Defining a holistic pain-relieving approach to wound care via a drug-free polymeric membrane dressing. *J Wound Care* 2011; 20(5): 250-4
10. Hayden JK, Cole BJ. The effectiveness of a pain wrap compared to a standard dressing on the reduction of post-operative morbidity following routine knee arthroscopy: A prospective randomized single-blind study. *Orthopedics* 2003; 26: 59-63
11. Sessions RC. Can a drug-free dressing decrease inflammation and wound pain? What does the evidence say? Poster IR-09. SAWC, September 2009.
12. Denyer J. Managing pain in children with epidermolysis bullosa. *Nurs Times* 2012; 108(29): 21-23
13. Scott A. Polymeric membrane dressings for radiotherapy-induced skin damage. *Br J Nurs (Oncology Suppl)* 2014; 23(10): S24-31
14. Wounds International. Optimising wellbeing in patients living with a wound. An international consensus. *Wounds International*, 2013
15. Benskin L. PolyMem WIC Silver Rope: a multifunctional dressing for decreasing pain, swelling, and inflammation. *Adv Wound Care* 2012; doi 10.1089/wound.2011.0285
16. Wilson D. New PolyMem Wic Silver Rope cavity filler dressing solves difficult tunnelling wound problem. Poster. SAWC, October 2008.
17. Benskin L. Extensive tunneling lower leg wounds with exposed tendons closed quickly using various PolyMem dressings. Poster. SAWC, October 2008.
18. Agathangelou C. An easy, effective and almost painless way to debride and heal arterial ulcers. Poster. EWMA, 2011.
19. Agathangelou C. Three years' experience of treating 46 painful diabetic foot ulcers with polymeric membrane dressings. Poster. EWMA, 2012
20. Rahman S, Shokri A. Total knee arthroplasty infections eliminated and rehabilitation improved using polymeric membrane dressing circumferential wrap technique: 120 patients at 12-month follow up. Poster. EWMA, 2013
21. Haik J, Weissman O, Demetris S, et al. Polymeric membrane dressings for skin graft donor sites. Six years' experience on 1200 cases. Poster. WUWUS, Japan, 2012
22. Tamir J, Haik J. Polymeric membrane dressings for skin graft donor sites: 4 years' experience on 800 cases. Poster. SAWC, October 2008
23. Skrinjar E, Duschek N, Bayer G et al. Efficiency of local wound treatment by combining polymeric membrane dressings and negative pressure wound therapy. Poster. EWMA, 2011
24. Vanwalleghem G. Four years' experience of effectively treating compartment syndrome by combining NPWT with polymeric interface layer. Poster. EWMA, 2013
25. Irving V, Bethell E, Burton F. Neonatal wound care: Minimising trauma and pain. *Wounds UK* 2006; 2(1): 33-41

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